
DOC022.97.93082

# TitraLab® KF1000 series workstations

10/2021, Edition 11

**Basic User Manual**  
**Manuel d'utilisation simplifié**  
**Manual básico del usuario**  
**Manual básico do usuário**  
基本用户手册  
基本取扱説明書  
기본 사용 설명서  
دليل المستخدم الأساسي

## Table of Contents

---

English.....	3
Français.....	28
Español.....	54
Português.....	81
中文.....	107
日本語.....	130
한국어.....	156
العربية.....	180

## Table of Contents

- |   |  |
|---|--|
| 1 <a href="#">Online user manual</a> on page 3  | 5 <a href="#">User interface and navigation</a> on page 17 |
| 2 <a href="#">Specifications</a> on page 3      | 6 <a href="#">Startup</a> on page 19                       |
| 3 <a href="#">General information</a> on page 3 | 7 <a href="#">Standard operations</a> on page 20           |
| 4 <a href="#">Installation</a> on page 8        | 8 <a href="#">Maintenance</a> on page 26                   |

## Section 1 Online user manual

This Basic User Manual contains less information than the User Manual, which is available on the manufacturer's website.

## Section 2 Specifications

Specifications are subject to change without notice.

Specification	Details
Dimensions (W x D x H)	22 x 40 x 36 cm (8.7 x 15.7 x 14.2 in.)
Weight	4 kg (8.8 lb)
Power requirements	100–240 VAC, 50/60 Hz
Main supply voltage fluctuation	±10% of nominal voltage
Altitude	2,000 m (6,562 ft) maximum
Operating temperature	15 to 35 °C (59 to 95 °F)
Relative humidity	20 to 80%, non-condensing
Storage temperature	–5 to 40 °C (23 to 104 °F)
Installation category	II
Pollution degree	2
Certifications	Safety IEC/EN 61010-1; EMC IEC/EN 61326-1
EMC requirements	This product is intended to be used in a domestic or basic electromagnetic environment
Warranty	1 year (EU: 2 years)

## Section 3 General information

In no event will the manufacturer be liable for direct, indirect, special, incidental or consequential damages resulting from any defect or omission in this manual. The manufacturer reserves the right to make changes in this manual and the products it describes at any time, without notice or obligation. Revised editions are found on the manufacturer's website.

### 3.1 Safety information

The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

Make sure that the protection provided by this equipment is not impaired. Do not use or install this equipment in any manner other than that specified in this manual.

### 3.1.1 Use of hazard information

<b>▲ DANGER</b>
Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.

<b>▲ WARNING</b>
Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.

<b>▲ CAUTION</b>
Indicates a potentially hazardous situation that may result in minor or moderate injury.

<b>NOTICE</b>
Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

### 3.1.2 Precautionary labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. A symbol on the instrument is referenced in the manual with a precautionary statement.

	This symbol, if noted on the instrument, references the instruction manual for operation and/or safety information.
	This symbol indicates that a risk of electrical shock and/or electrocution exists.
	This symbol indicates the presence of devices sensitive to Electro-static Discharge (ESD) and indicates that care must be taken to prevent damage with the equipment.
	Electrical equipment marked with this symbol may not be disposed of in European domestic or public disposal systems. Return old or end-of-life equipment to the manufacturer for disposal at no charge to the user.

### 3.1.3 Compliance and certification

<b>▲ CAUTION</b>
This equipment is not intended for use in residential environments and may not provide adequate protection to radio reception in such environments.

#### **Canadian Radio Interference-Causing Equipment Regulation, ICES-003, Class A:**

Supporting test records reside with the manufacturer.

This Class A digital apparatus meets all requirements of the Canadian Interference-Causing Equipment Regulations.

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### **FCC Part 15, Class "A" Limits**

Supporting test records reside with the manufacturer. The device complies with Part 15 of the FCC Rules. Operation is subject to the following conditions:



1. The equipment may not cause harmful interference.
2. The equipment must accept any interference received, including interference that may cause undesired operation.

Changes or modifications to this equipment not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference, in which case the user will be required to correct the interference at their expense. The following techniques can be used to reduce interference problems:

1. Disconnect the equipment from its power source to verify that it is or is not the source of the interference.
2. If the equipment is connected to the same outlet as the device experiencing interference, connect the equipment to a different outlet.
3. Move the equipment away from the device receiving the interference.
4. Reposition the receiving antenna for the device receiving the interference.
5. Try combinations of the above.

## **3.2 Product overview**

The instrument operates with digital and analog sensors. Measurement applications are installed on the instrument to automate the measurement process. Instructions show on the display when user intervention is required.

1 Keypad	6 Syringe protection cover	11 Pump 2 input/output (waste)
2 Display	7 Syringe input/output	12 Pump 1 input/output (solvent)
3 USB port <sup>1</sup>	8 Tube clips	13 Pump access cover
4 Sample stopper	9 Syringe	
5 Beaker	10 Sensor holder	


**Table 1 Instrument configuration**

Model	Syringes	Pumps
KF1121	1	2

<sup>1</sup> A second USB port is found on the rear of the instrument, but the instrument recognizes only one USB storage device connection at a time.

### 3.3 Instrument connections

Use the USB port on the side of the instrument for the USB applications key supplied with the instrument. Use the USB port on the rear of the instrument to connect to a printer, mouse, keyboard or a USB hub.


1 24 V external power supply port	4 Not used	7 USB port
2 Sensor 1 port	5 Not used	8 Ethernet port
3 Sensor 2 port	6 Serial port	

### 3.4 Product components

Make sure that all components have been received. Refer to the packing list in the box. If any items are missing or damaged, contact the manufacturer or a sales representative immediately.


**Figure 1 Contents of the instrument box**

*Note: For identification, the outlet tube from the syringe is fitted with a blue marker.*

1 Instrument	3 Power cord
2 Sample stopper	4 Power supply

**Figure 2 Contents of the application box**

1 KF standard beaker	6 Syringe holding ring
2 KF white PTFE conical adapter	7 Syringe
3 Magnetic stir bars	8 Glass bottle
4 USB applications key	9 Bottle caps (3 x GL45)
5 Sensor	10 Desiccant cartridges (3x, filled with a molecular sieve)

## Section 4 Installation

### ⚠ CAUTION

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

### NOTICE

This is a class A product. There may be potential difficulties in ensuring electromagnetic compatibility in other environments, due to conducted as well as radiated disturbances. In a domestic environment this product may cause radio interference in which case the user may be required to take adequate measures.

### NOTICE

Network and access point security is the responsibility of the customer that uses the wireless instrument. The manufacturer will not be liable for any damages, inclusive however not limited to indirect, special, consequential or incidental damages, that have been caused by a gap in, or breach of network security.

### 4.1 Installation guidelines

- This instrument is for indoor use only.

- The power supply connector on the rear panel must be easily accessible so the power can be disconnected quickly in case of emergency.
- Keep the instrument away from temperature extremes, including heaters, direct sunlight and other heat sources.
- Put the instrument on a stable and level surface in a well ventilated place.
- Make sure that there is at least 15 cm (6 in.) of space on all sides of the instrument to prevent electrical parts from overheating.
- Do not operate or keep the instrument in dusty, damp or wet locations.
- Always keep the surface of the instrument and all accessories dry and clean.

## 4.2 Connect to AC power

<b>⚠ DANGER</b>	
	Electrocution hazard. If this equipment is used outdoors or in potentially wet locations, a Ground Fault Circuit Interrupt (GFCI/GFI) device must be used for connecting the equipment to its main power source.
<b>⚠ CAUTION</b>	
	Electrical shock and fire hazards. Make sure that the supplied cord and non-locking plug meet the applicable country code requirements.
<b>⚠ WARNING</b>	
	Fire hazard. Use only the external power supply that is specified for this instrument.

1. Connect the power cord to the power supply.
2. Connect the power supply to the instrument (refer to [Instrument connections](#) on page 7).
3. Connect the power cord to an electrical outlet.

## 4.3 Install the syringe







Before syringe installation, set the instrument power to on. Push the power button on the front of the instrument. Make sure that the startup sequence shows on the display. The syringe holder lowers to its operating position.



*Note: Ignore any warning messages related to missing applications that show on the display.*

The sensor holder has two positions: one over the magnetic stirrer and the second at 180° to the right. Move the sensor holder away from the instrument to the second position.

Refer to the illustrated steps that follow.


<b>⚠ CAUTION</b>	
In step 6, tighten the syringe using the metallic part at the top. Do not hold the glass section of the syringe. Do not tighten too much.	


#### 4.4 Install the sensor holder accessories

Add the stir bar to the beaker, and then attach the beaker to the sensor holder. Install the desiccant cartridge and the sample stopper. Refer to the illustrated steps that follow.

#### 4.5 Prepare the tubes

Remove any bends in the end of the tubes.


## 4.6 Connect the tubes

Arrow symbols identify the inlet and outlet ports for the syringe and the pump connections. The “up” arrow is the outlet port. The “down” arrow is the inlet port. Turn the tube connectors on the inlet and outlet ports of the syringe and pump until they click.

The suction tube that drains into the waste bottle must be installed on the bottom of the beaker for good drainage.

### NOTICE

The diffusion tip and tube holder on the outlet tube, are pre-installed at the optimal positions. Do not change the position of the diffusion tip or the tube holder.


## 4.7 Install the sensor

### 4.7.1 Connect the sensor

Use a conical adapter to hold the sensor tightly in the sensor holder. Make sure the conical adapter is installed correctly. This makes sure of a secure water-tight fitting in the sensor holder and the tip of the sensor will be in the correct position in the beaker.

Connect the sensor to an available sensor port on the rear of the instrument. After the sensor is connected, make sure that the sensor icon shows in the banner at the top of the display.


Refer to the illustrated steps that follow.

#### 4.7.2 Check the titrant tube and sensor installation

The titrant must be added where there is most turbulence and as far away from the sensor to allow the reagent time to react. Check that the titrant tube and the sensor are positioned correctly. Refer to [Figure 3](#), where the arrows show the direction of the sample flow.

**Figure 3** Titrant tube and sensor positions

1 Titrant tube	2 Sensor	3 Magnetic stir bar
----------------	----------	---------------------

#### 4.8 Install the titrant and the reagent

##### ▲ CAUTION

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

##### ▲ CAUTION




Chemical exposure hazard. Dispose of chemicals and wastes in accordance with local, regional and national regulations.

Put a filled desiccant cartridge into the adapter on the titrant and solvent bottle caps.


**Note:** The molecular sieve in the filled cartridges must be regenerated regularly depending on the humidity of the laboratory. Refer to [Regenerate the molecular sieve](#) on page 27.

Loosen the tube connector on the bottle cap. Push the inlet tube through the connector. Make sure that the end of the tube is at the bottom of the bottle. Tighten the connector on the bottle cap.

Pump 1 is used to fill the measurement cell with solvent. Pump 2 is used to empty the measurement cell. Refer to the illustrated steps that follow.

**Figure 4 Bottles identification**


1 Titrant

2 Waste

3 Solvent

## 4.9 Tidy the work area


Attach the tubes to the instrument with the clips on the electrovalve and the sensor holder. Refer to the illustrated steps that follow.

## Section 5 User interface and navigation


### 5.1 Keypad

1 Power	3 Selection keys	5 Home key
2 Printer	4 Navigation keys	

Key	Description
Power	Sets the instrument power to on or off. Push the key for 2 seconds to set the power to off.
Printer	The printer key only operates if a printer is connected to the instrument. Sends the data currently shown on the display to a connected printer when pushed. A sound is heard if the current display cannot be printed. A graph is automatically printed at the end of the measurement, if the option is selected ( <b>Settings &gt; Options</b> ).
Selection keys (contextual)	Used to select options shown above them in the footer bar. The options available are applicable to the current operation (e.g. calibration, measurement, etc.).
Navigation keys	Scrolls through menus and data, enter numbers and letters, enter checkbox settings and set options for the syringe and the pump.
Home	Push this key at any time to go directly to the home screen. A sound is heard if the key is disabled (e.g. during a calibration or measurement).

## 5.2 Home screen

1 Available options from this screen	5 Screen name
2 A second row of available options from this screen	6 Information icons (refer to <a href="#">Table 2</a> )
3 KF applications	7 Option available by pushing the selection key below
4 Screen icon	8 Arrow keys available for use in the screen

[Table 2](#) shows the information icons that can be shown in the header bar.

**Table 2 Information icons**

Icon	Description
	A printer is connected to the instrument
	A USB key is connected to the instrument
	A sensor is connected to sensor port 1
	A sensor is connected to sensor port 2
	The data log file is full. Refer to <a href="#">Manage the data log</a> on page 25 for the options available to manage the data log file.

## Section 6 Startup

### ▲ CAUTION

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

### ▲ CAUTION

Chemical exposure hazard. Dispose of chemicals and wastes in accordance with local, regional and national regulations.

### ▲ CAUTION

Personal injury hazard. Never use the instrument without the syringe cover installed.

## 6.1 Configure the instrument

1. From the main menu, select **Settings**.
2. Select an option, then push **Select**.

Option	Description
<b>Applications</b>	Changes, makes copies, exports and removes application data. Make sure that the duplication function does not make more than five applications.
<b>Operators</b>	Adds, changes and removes operators.
<b>Date + Time</b>	Sets the instrument date and time.
<b>Brightness</b>	Sets the brightness of the display.
<b>Sounds</b>	Sets the sound options.
<b>Language</b>	Sets the language.
<b>Network</b>	Give a name to the instrument. If a printer is connected, this name is printed on the hard-copy output. Restart the instrument if the name is changed.
<b>Info</b>	Shows information about the instrument and the attached hardware.
<b>Restore Defaults</b>	Sets the instrument to the default configuration.
<b>Options</b>	Sets the application parameters view to expert mode. When the instrument is set to off, sets the syringe to empty into the titrant bottle. Changes the temperature display unit. Set the cell standby parameter (refer to <a href="#">Cell standby</a> on page 19). Prints the measurement and derivative curves if a printer is connected. Specify if a balance is connected.
<b>Security</b>	Change the password and specify which options are password protected.

3. Push **Back**.

### 6.1.1 Cell standby

Make sure that the measurement cell is prepared before selecting an application, because the cell standby starts automatically.

Cell standby dries the measurement cell automatically using a regulated addition of the titrant as defined in the application. When the cell standby option is active:

- Cell standby starts automatically as soon as an application is selected from the home screen
- The stirring speed is the speed specified for the application

- Cell standby information is shown at the bottom of the application start screen
- After the weight of the sample is validated at the end of the measurement procedure, the cell conditioning automatically starts again
- Cell standby continues until the next sample is selected and a new measurement procedure starts
- If the **Cell management** option is selected at the end of the measurement procedure, cell standby is paused and starts again when the cell management procedure completes
- The measurement procedure starts when **Start** is selected. If the cell is dry, cell conditioning is not needed and the user is asked to add the sample. If the cell is not dry, cell conditioning starts

## 6.2 Install the applications

Use the supplied USB key to install the applications. The instrument can install a maximum of five applications. The installed applications are shown on the top line of the display.

1. Push **Home** to go to the main menu.
2. Connect the USB key to the USB port on the side of the instrument. The applications on the USB key show on the display.
3. Push the arrow keys to highlight and select an application to install. Push the left or right arrow key to select it. Do this step again to select additional applications to install.
4. Push **Import** to install the selected applications.
5. Push **OK** to complete the installation. The installed applications show on the home screen.

*Note: To install more applications, push **Home** to go to the home screen, then remove the USB key and reconnect it.*

## 6.3 Prepare the instrument for measurement

1. From the home screen, select **Purge**. All attached devices are listed.
2. Select **All elements** to purge all the attached devices, or select one device to purge. Push **Select**. Air is removed from the device and filled with liquid from the bottle.
3. Push **OK** when the operation has completed.
4. Make sure that there are no air bubbles in the device. Do step 2 again if there are any air bubbles.
5. Select the next device to purge if individual devices are being selected.
6. Push **Exit** when all the tubes are filled with reagent and the device has no air bubbles.

*Note: If a few small air bubbles can be seen on the inner wall and/or piston of the syringe, they can be left without affecting system performance.*

## Section 7 Standard operations

### ▲ CAUTION

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

### ▲ CAUTION

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

### ▲ CAUTION

Chemical exposure hazard. Dispose of chemicals and wastes in accordance with local, regional and national regulations.

## ▲ CAUTION

Personal injury hazard. Never use the instrument without the syringe protection cover in place.

## ▲ CAUTION

Chemical exposure hazard. Never remove the stir bar from the beaker before the end of a titration.

## 7.1 Calibration

### 7.1.1 Calibrate the titrant

1. Make sure that the glass beaker, the sensor holder and all related parts are clean and dry. Put a magnetic stir bar into the beaker.
2. Attach the beaker to the sensor holder. Make sure that the beaker fits tightly into the sensor holder.
3. Read the related "Application Note" from the USB applications key for more instructions.
4. Fill the beaker with the applicable amount of solvent. Refer to the application note and [Cell management](#) on page 26 to adapt the volume where applicable. Make sure that there is no liquid spill. Make sure that the initial level is enough so that the sensor is correctly installed in the sample. Do not put too much solvent in the beaker.
5. From the home screen, select **Calibration**.
6. If more than one application contains a titrant calibration method, push the up and down arrow keys to highlight the application to use, then push **Select**.
7. Calibration information shows on the display. If necessary, select an icon for more information or to change some data.

Option	Description
<b>Information</b>	Shows more information about the calibration.
<b>Operator</b>	Changes the operator ID. Select from a list of applicable operators.

8. Make sure that the icon at the bottom of the display  is highlighted. Do the instructions that show on the display adjacent to this icon. Make sure that the tubes and sensor are correctly aligned.
9. Push **Start** to start the calibration. Calibration data shows on the display.
10. If the default stirring speed needs to be adjusted, push the up and down arrow keys to increase or decrease the speed.  
*Note: This adjustment only applies to the current operation. The standard default stirring speed for the application is not changed.*
11. The cell drift threshold, stability bar, sample temperature and measurement duration are shown at the bottom of the display. The stability bar shows the stability of the signal and ranges in color from red (unstable), through yellow, to green (stable). The cell drift value is shown in the middle of the display. These messages can be displayed under the cell drift value:

Message	Description
<b>Excess of water</b>	The measurement cell contains too much water. Titrant continues to be added until the cell is dry.
<b>Excess of titrant</b>	The measurement cell contains too much titrant. No more titrant is added.
<b>Risk of cell overflow</b>	The measurement cell is full and there is a risk of liquid spill. The procedure stops.

*Note: The stability bar or these messages are displayed at the bottom of the display during the calibration procedure.*

12. Two options are available during the cell conditioning procedure:

Option	Description
--------	-------------

<b>Stop</b>	Aborts the measurement and no results are calculated. If selected during the <b>Replicate Sample</b> option, then all data in the series is lost.
-------------	---

<b>Skip</b>	Stops the cell conditioning and goes directly to the next step in the procedure. Results are then calculated from measurement data available before <b>Skip</b> is selected. Results can be less accurate using this option.
-------------	--

13. Weigh the standard in its container and write down the result.

14. When the cell conditioning is complete add the standard. Remove the sample stopper and put the correct amount of standard into the beaker with an applicable utensil. Make sure that all of the standard is put into the center of the beaker and that none of it spills onto the walls of the beaker.

15. Put the sample stopper back into the sensor holder. Make sure that the sample stopper fits tightly into the sensor holder.

16. Weigh the empty sample container and subtract this value from the weight of the standard in its container (refer to step 13). The result is the exact weight of the sample analyzed. Write down this value.

17. Push **Start** to start the moisture content analysis. Calibration data shows on the display.

*Note: The moisture content analysis will start automatically after the addition of sample if the **Injection Autodetect** parameter is set to **Yes**.*

18. If the default stirring speed needs to be adjusted, push the up and down arrow keys to increase or decrease the speed.

*Note: Changing the default stirring speed can effect the measurement result.*

19. Two options are available during the calibration procedure:

Option	Description
--------	-------------

<b>Stop</b>	Aborts the calibration and no results are calculated. If selected during the <b>Replicate Sample</b> option, then all data in the series is lost.
-------------	---

<b>Skip</b>	Stops the calibration and goes directly to the next step in the procedure. Results are then calculated from calibration data available before <b>Skip</b> is selected. Results can be less accurate using this option.
-------------	--

20. When the calibration is complete, enter the exact weight of the sample analyzed (refer to step 16). If necessary, push the arrow keys to change the value shown on screen. Push **Select** to accept the value.

21. When the calibration is complete, push the left and right arrow keys to see the different calibration views.

22. Push **Reject** or **Continue**.

Option	Description
--------	-------------

<b>Reject</b>	Rejects the calibration. Select <b>Cancel</b> to go back to the result display or <b>Confirm</b> to reject the calibration. If this is the first calibration, select <b>Confirm</b> to reject the calibration and use the default or previous calibration values. If this is a <b>Replicate Standard</b> calibration, select <b>Confirm</b> to reject only the current calibration in the series.
---------------	---

<b>Continue</b>	Select from one of these options:
-----------------	-----------------------------------

- **Replicate Standard:** Do the calibration again using the same standard
- **Save & Exit:** Keep the calibration results and exit the calibration procedure
- **Reject & Exit:** Reject the calibration results and use the default or previous calibration values, and exit the calibration procedure

## 7.2 Get a sample measurement

Use this option to get sample measurements with one of the installed applications.

1. Make sure that the glass beaker, the sensor holder and all related parts are clean and dry. Put a magnetic stir bar into the beaker.
2. Attach the beaker to the sensor holder. Make sure that the beaker fits tightly into the sensor holder.
3. Read the related "Application Note" from the USB applications key for more instructions.
4. Fill the beaker with the applicable amount of solvent. Refer to the application note and [Cell management](#) on page 26 to adapt the volume where applicable. Make sure that there is no liquid spill. Make sure that the initial level is enough so that the sensor is correctly installed in the sample. Do not put too much solvent in the beaker.
5. From the home screen, select the measurement application, then push **Select**.
6. If the **Cell standby** option has not been selected (refer to [Cell standby](#) on page 19), continue at step 8.
7. The cell standby operation starts. The stirring speed is that specified for the application and cannot be changed. The stability bar and the cell drift value are shown at the bottom of the display. The stability bar shows the stability of the signal and ranges in color from red (unstable), through yellow, to green (stable). The cell drift value can be replaced by these messages:

Message	Description
<b>Excess of water</b>	The measurement cell contains too much water. Titrant continues to be added until the cell is dry.
<b>Excess of titrant</b>	The measurement cell contains too much titrant. No more titrant is added.
<b>Risk of cell overflow</b>	The measurement cell is full and there is a risk of liquid spill. The procedure stops.

*Note: These messages can be displayed under the cell drift value in the middle of the display during the cell conditioning procedure. The stability bar or these messages are displayed at the bottom of the display during the measurement procedure.*

8. Application information shows on the display. If necessary, select an icon for more information or to change some data.

Option	Description
<b>Information</b>	Shows more information about the application.
<b>Operator</b>	Changes the operator ID. Select from a list of applicable operators.
<b>Sample</b>	Sample Name: Change the specified name of the sample. Type: Push the left and right arrow keys and select the sample type ( <b>Sample</b> , <b>QC</b> or <b>Define blank</b> ) to be used for the measurement. If <b>Define blank</b> has been selected before, two more sample types are available ( <b>QC with blank</b> and <b>Sample with blank</b> ).

9. Make sure that the icon at the bottom of the display  is highlighted. Do the instructions that show on the display adjacent to this icon. Make sure that the tubes and sensor are correctly aligned.
10. Push **Start** to start the measurement. Measurement data shows on the display.
11. If the **Cell standby** option has been selected (refer to [Cell standby](#) on page 19), and the measurement cell is dry, cell conditioning is not needed and the procedure continues at step 15.

12. If the default stirring speed needs to be adjusted, push the up and down arrow keys to increase or decrease the speed.

*Note: Changing the default stirring speed can effect the measurement result. Any adjustment only applies to the current operation. The standard default stirring speed for the application is not changed.*

13. The cell drift threshold, stability bar, sample temperature and measurement duration are shown at the bottom of the display. The stability bar shows the stability of the signal and ranges in color from red (unstable), through yellow, to green (stable). The cell drift value is shown in the middle of the display. These messages can be displayed under the cell drift value:

Message	Description
<b>Excess of water</b>	The measurement cell contains too much water. Titrant continues to be added until the cell is dry.
<b>Excess of titrant</b>	The measurement cell contains too much titrant. No more titrant is added.
<b>Risk of cell overflow</b>	The measurement cell is full and there is a risk of liquid spill. The procedure stops.

*Note: The stability bar or these messages are displayed at the bottom of the display during the measurement procedure.*

14. Two options are available during the cell conditioning procedure:

Option	Description
<b>Stop</b>	Aborts the measurement and no results are calculated. If selected during the <b>Replicate Sample</b> option, then all data in the series is lost.
<b>Skip</b>	Stops the cell conditioning and goes directly to the next step in the procedure. Results are then calculated from measurement data available before <b>Skip</b> is selected. Results can be less accurate using this option.

15. Weigh the sample in its container and write down the result.

16. When the cell conditioning is complete add the sample. Remove the sample stopper and put the correct amount of sample into the beaker with an applicable utensil. Make sure that all of the sample is put into the center of the beaker and that none of it spills onto the walls of the beaker.

17. Put the sample stopper back into the sensor holder. Make sure that the sample stopper fits tightly into the sensor holder.

18. Weigh the empty sample container and subtract this value from the weight of the sample in its container (refer to step 15). The result is the exact weight of the sample analyzed. Write down this value.

19. Push **Start** to start the moisture content analysis. Measurement data shows on the display.

*Note: The moisture content analysis will start automatically after the addition of sample if the **Injection Autodetect** parameter is set to **Yes**.*

20. If the default stirring speed needs to be adjusted, push the up and down arrow keys to increase or decrease the speed.

*Note: Changing the default stirring speed can effect the measurement result.*

21. Two options are available during the measurement procedure:

Option	Description
<b>Stop</b>	Aborts the measurement and no results are calculated. If selected during the <b>Replicate Sample</b> option, then all data in the series is lost.
<b>Skip</b>	Stops the measurement and goes directly to the next step in the procedure. Results are then calculated from measurement data available before <b>Skip</b> is selected. Results can be less accurate using this option.

22. When the measurement is complete, enter the exact weight of the sample analyzed (refer to step 18). If necessary, push the arrow keys to change the value shown on screen. Push **Select** to accept the value.
23. If the **Cell standby** option has been selected (refer to **Cell standby** on page 19), the standby operation starts again automatically.
24. Push the arrow keys to see the different measurement views or push **Next** for the options that follow:

Option	Description
<b>Replicate sample</b>	Use this option to start the same titration on the same sample. This is used to study the repeatability by successively analyzing several parts of the same sample. At the end of each measurement, a window shows the average value, the standard deviation and the relative standard deviation.
<b>New sample</b>	Use this option to start the same titration on a new sample. No standard deviation or relative standard deviation measurements will be done.
<b>Cell management</b>	Use this option to fill or empty the measurement cell using the pumps. Refer to <b>Cell management</b> on page 26.

25. Push **Exit** to go back to the home screen.

### 7.3 Syringe verification

Use the syringe verification option to examine the syringe accuracy based on the ISO 8655 standard. Make sure to fully rinse the syringe with *methanol dry* before the test to remove titrant. Make sure to fully rinse the syringe with *methanol dry* after the test to remove water.

*Note: A license key is necessary to enable the syringe verification option on the instrument. Contact the manufacturer or a sales representative to get a license key.*

1. From the home screen, select Settings, then push **Select**.
2. From the settings screen, select Options, then push **Select**.
3. Push the left or right arrow key to select the expert mode, then push **Done**.
4. From the home screen, select Maintenance, then push **Select**.
5. From the maintenance screen, select Syringe verification, then push **Select**.
6. When the option is selected for the first time, enter the license key.
7. Follow the procedure supplied with the syringe verification kit. For additional information, refer to the expanded version of this manual.

### 7.4 Manage the data log

To select data to view, delete or export, specify data filters

1. From the home screen, select **Data log**.
2. Select an option, then push **Select**.

Option	Description
<b>View data log</b>	Views measurement data. Select individual lines of data to view more content.
<b>Export data log</b>	Exports measurement data from the system to an external device. Preview data selection before it is exported. Make sure that an external device is connected to the instrument (e.g., a USB key, external hard drive, etc.).
<b>Delete data log</b>	Removes measurement data from the system. Previews data selection before it is removed.

- Specify the data selection parameters. Push the left and right arrow keys to make a selection. Push the up and down arrow keys to select an option.

Option	Description
<b>Result type</b>	Sets the type of result available.
<b>Application</b>	Sets the available applications.
<b>Date</b>	Sets the date range.
<b>Operator</b>	Sets the available operators.

- If **View data log** was the selected option, push **View** to see the selected data.
  - Push the up and down arrow keys to select a line of data and push **Detail** to see more data
  - If only one **Application** is selected push the up and down arrow keys to select a line of data and push **Detail** to see more data, or push the left and right arrow keys to see related graphs
- If **Export data log** or **Delete data log** was the selected option, push **Preview** to see the selected data, then push **Export** or **Delete** to start the procedure.

## 7.5 Purge

Use this procedure to remove air bubbles from the system. Refer to [Prepare the instrument for measurement](#) on page 20 for instructions.

## 7.6 Cell management

Use this option to fill or empty the measurement cell with the pumps. These procedures stop automatically after a pre-determined time. The time remaining is shown against **Security time** on the display.

- From the home screen, select **Cell management**.
- Push the left and right arrow keys to make a selection.

Option	Description
<b>Fill cell (pump 1)</b>	Use this option to fill the measurement cell with reagent from pump 1.
<b>Empty cell (pump 2)</b>	<div style="text-align: center; background-color: #FFD700; padding: 5px;"><b>⚠ CAUTION</b></div> <p>Chemical exposure hazard. Before selecting this option, examine the level of liquid in the waste bottle. Make sure that there is enough space in the bottle for all of the liquid from the measurement cell.</p> <p>Use this option to empty the measurement cell with pump 2.</p>
<b>Stop pumps</b>	Use this option to stop the <b>Fill cell</b> and <b>Empty cell</b> procedures.
<b>Stir</b>	Use this option to start or stop the magnetic stirrer. This option can also be used when the <b>Fill cell</b> and <b>Empty cell</b> procedures are operating.

- Push **Exit** to go back to the home screen.

## Section 8 Maintenance

<b>⚠ CAUTION</b>	
	Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

## **NOTICE**

Do not disassemble the instrument for maintenance. If the internal components must be cleaned or repaired, contact the manufacturer.

### **8.1 Clean the instrument**

## **NOTICE**

Never use flammable or corrosive solvents to clean any part of the instrument. Use of these solvents can degrade the environmental protection of the instrument and may void the warranty.

Clean the exterior surface with a moist cloth or with a mixture of water and mild detergent. Dry with a soft cloth.

### **8.2 Clean the sensor**

Refer to the documentation delivered with the sensor.

### **8.3 Regenerate the molecular sieve**

It is recommended to regenerate the molecular sieve weekly, but this will vary depending on the humidity of the laboratory. A humidity indicator, such as color changing silica gel crystal, can be used to show when regeneration is necessary.

To regenerate the molecular sieve, dry it in an oven at a temperature of 300 °C (572 °F) for at least 4 hours. When cool, it can be kept in a glass bottle with an airtight seal.

### **8.4 Maintenance menu**

Refer to the Full User Manual, which is available as a download from our website.

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P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
(800) 227-4224 (U.S.A. only)