
DOC012.97.90328

# SIP 10

05/2011, Edition 2

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## Specifications

Subject to change without notice.

Performance specifications	
Cell path length	Dual pathlength Pour-Thru Cell, plastic, 1 inch/ 1 cm
Wavelength range	340–900 nm
Rinsing volume	Minimum of 25 mL
Flow speed	1 mL/sec
Interface	USB type A
Power supply	Via USB cable, 530 mA, 5 V
Ambient storage requirements	–10–60 °C (14–140 °F), 85% relative humidity (non-condensing)
Ambient operating requirements	10–40 °C (50–104 °F), 80% relative humidity (non- condensing)
Dimensions	120 x 85 x 200 mm
Weight	0.5 kg
Enclosure rating	IP30

## General Information

### Safety Information

Please read the entire user manual carefully before the device is unpacked, set up or put into operation. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

To make sure that the safeguards provided by this device are not impaired, do not use or install this device in any manner other than that specified in this user manual.

### Hazard Information in this User Manual

#### **DANGER**

Indicates a potentially or imminently hazardous situation that, if not avoided, results in death or serious injury.

#### **WARNING**

Indicates a potentially or imminently hazardous situation that, if not avoided, could result in death or serious injury.

#### **CAUTION**

Indicates a potentially hazardous situation that may result in minor or moderate injury.

#### **NOTICE**

Indicates a situation that, if not avoided, could result in damage to the instrument. Information that requires special emphasis.

**Note:** Information that supplements points in the main text.

## Precautionary Labels

Read all labels and tags attached to the device. Personal injury or damage to the device could occur if this information is not observed.

	This symbol is a warning triangle. Follow all safety instructions preceded by this symbol to prevent potential injury. If this symbol is attached to the device, it refers to information in the operation and/or safety sections of the user manual.
	This symbol may be found on an enclosure or barrier within the product and indicates a risk of electric shock and/or death by electrocution.
	Electrical devices marked with this symbol may not be disposed of via European domestic or public disposal systems after 12 August 2005. In accordance with current regulations (EU Directive 2002/96/EC), users within the EU must now return old electrical devices to the manufacturer for disposal. The user bears no costs for this process. <b>Note:</b> <i>Instructions on the correct disposal of all (marked and unmarked) electrical products supplied or manufactured by HACH LANGE can be obtained from your local HACH LANGE sales office.</i>

## Chemical and Biological Safety

### **WARNING**

Chemicals can cause serious injuries if handled without the relevant protection. Always wear personal protective equipment to detach and replace tubes in accordance with the safety data sheets applicable to these chemicals.

### **CAUTION**

Chemical hazard. Do not use the Pour-Thru Cell in tests that require the use of organic solvents such as alcohols, toluene, chloroform, trichloroethane or cyclohexanone. These solvents may not be compatible with the plastic components of the Pour-Thru Cell, which creates the potential for damage to the device and chemical exposure for the analyst.

Normal operation of the device may involve the use of hazardous chemicals or biologically harmful samples.

- Observe all cautionary information printed on the original solution containers and safety data sheet prior to their use.
- Dispose of all consumed solutions in accordance with national regulations and laws.
- Select the type of protective equipment suitable to the concentration and quantity of the dangerous material at the respective work place.
- Always rinse the cell and tubes thoroughly after use with deionized water.

## Product Overview

The SIP 10 Sipper module is an optional accessory for the DR 3900 spectrophotometer. The SIP 10 serves to improve measurement accuracy. The same optical characteristics apply to both the zero compensation and the measurements. All test solutions flow through a cell and any errors due to optical differences between redundant cells are eliminated. A peristaltic pump pumps a constant amount of fluid through the Pour-Thru Cell.

## Scope of Delivery

Please check that the order is complete. If anything is missing or damaged, please contact the manufacturer or distributor immediately.

- SIP 10 Sipper module, pump tubing pre-assembled
- Dual pathlength Pour-Thru Cell, 1 inch/1 cm
- USB connection cable
- Tygon® tubing, 1.70 m (5.6 ft)
- Stainless steel tube adapter
- User manual, multilingual

## Installation

### NOTICE

Make sure that there are no kinks in the various tubes.  
Check the fill level of the waste vessel on a regular basis.

### NOTICE


There may always be residual fluid in the tubes. Always wear suitable protective clothing and goggles to change the tubes or perform maintenance work in order to prevent contact with the skin or eyes.

### NOTICE

Connect the device directly to the photometer via the USB cable.  
Do not use a USB hub connection.

1. Position the SIP 10 to the left of the photometer.
2. Clean the Pour-Thru Cell with a lint-free cloth.
3. Insert the Pour-Thru Cell into the cell compartment of the photometer in the direction indicated by the arrow.
4. Cut a piece of tubing no longer than 40 cm (1.3 ft) to serve as the feed tube.
5. Connect one end of the feed tube to the IN connection of the Pour-Thru Cell.
6. Connect the other end of the feed tube to the stainless steel tube adapter and place this in the zero-solution/sample vessel.
7. Cut a piece of tubing approx. 30 cm (1 ft) long to serve as the connection tube.
8. Connect the connection tube to the OUT connection of the Pour-Thru Cell and the input connection of the peristaltic pump.
9. Use the remaining tubing as a discharge tube, and connect one end to the output connection of the peristaltic pump.
10. Place the other end of the discharge tube in a suitable waste vessel.
11. Connect the SIP 10 USB cable to a USB port on the photometer.

Figure 1 Fully installed SIP 10 module


1 SIP 10 Sipper module	6 Pour-Thru Cell
2 Discharge tube	7 Feed tube
3 Peristaltic pump	8 Photometer
4 Connection tube	9 Stainless steel tube adapter
5 USB cable	

## Operation

### SIP 10 setup

When a SIP 10 module is connected to the photometer, an additional **SIPPER** button is shown in the toolbar.

1. Press **SIPPER** on the toolbar.  
A Sipper options menu opens.
2. Press **ON** to modify the settings for a Sipper cycle.

- **SIP TIME:** The sip time determines the sample volume pumped through the cell.  
To make sure that there is a complete exchange of fluid in the Pour-Thru Cell, the minimum time setting is 25 seconds. Default setting: 25 seconds. The flow speed is 1 mL/sec.
- **SETTLE TIME:** The settle time defines the duration of the standstill phase between the pump process and measurement process. This time period allows any air bubbles that have formed during the pump process to disappear and any sample turbulence to cease.  
Default setting: 5 seconds
- **PURGE TIME:** The purge time determines the volume of rinse agent that is pumped through the cell after each measurement. This additional purge cycle can be deactivated.  
Default setting: off

- PURGE START:** The start of the purge cycle can be triggered automatically after the measurement process or manually. To start the cycle manually, press **PURGE** to trigger the purge cycle.  
 Default setting: auto

**Note:** The **PURGE START** button is only active when the purge time setting is ON.

- Configure the required settings and press **OK** to confirm.

## Execution of a Stored Program


- Access the required test via **STORED PROGRAMS**.
- Place the feed tube in the zero solution and press **ZERO**.

The zero solution is sucked in over the selected sip time. Once the settle time has passed, the zero measurement is taken.

**Note:** The remaining time is shown on the display. Press **CANCEL** to abort the operation.

**Note:** If the **PURGE START: AUTO** setting has been activated in the Sipper options menu, the rinse agent is sucked in automatically once the measurement process is complete.

The zero measurement is shown on the display.

- Place the feed tube in the sample solution and press **READ**.


The sample solution is sucked in over the selected sip time. Once the settle time has passed, the measurement is taken.

**Note:** The remaining time is shown on the display. Press **CANCEL** to abort the operation.

**Note:** If the **PURGE START: AUTO** setting has been activated in the Sipper options menu, the rinse agent is sucked in automatically once the measurement process is complete.

**Note:** Either deionized water or the next sample can be used as a rinse agent.

The measurement is shown on the display.

For the tests listed in [Table 1](#), the Pour-Thru Cell requires a purge cycle with deionized water between samples.

**Table 1 Tests that require a purge cycle**

Aluminum, Aluminon	Chlorine Dioxide, LR	Cobalt, PAN
Copper, Porphyrin	Hardness, Calmagite	Manganese, LR, PAN
Nickel, PAN	Nitrate, MR	Nitrate, HR

The Pour-Thru Cell can also be used for the Nessler Method for ammonia and TKN if cleaned correctly. Pour a few sodium thiosulfate crystals into the cell to clean it. Rinse out the crystals with deionized water.

The Pour-Thru Cell cannot be used for the tests listed in [Table 2](#) due to potential problems caused by chemical incompatibility or other complications. Use the sample cell specified in the procedural instructions for these tests.

**Table 2 Tests for which the Pour-Thru Cell is not suitable**

Aluminum ECR	Arsenic	Barium	Boron, Carmine
Cyanuric Acid	Fluoride	Formaldehyde	Lead, LeadTrak
Mercury	Nickel, Heptoxime	Nitrite, HR	PCB
Phenols	Potassium	Selenium	Silver
Suspended Solids	Sulfate	TPH	Volatile Acids
Zinc	Surfactants, Anionic (Detergents)		

The Pour-Thru Cell can also be used for other tests. However, please note that the Pour-Thru Cell requires a minimum sample volume of 25 mL to make sure that the previous solution is fully rinsed out of the cell. Make sure that you increase the reagent and sample volumes accordingly.

## Execution of a Single or Multi-Wavelength Measurement


1. Press **SINGLE WAVELENGTH** or **MULTI-WAVELENGTH**.
2. Press **OPTIONS** to adjust the settings within single/multi-wavelength mode.
3. Place the feed tube in the zero solution and press **ZERO**.

The zero solution is sucked in over the selected sip time. Once the settle time has passed, the zero measurement is taken.

**Note:** The remaining time is shown on the display. Press **CANCEL** to abort the operation.

**Note:** If the **PURGE START: AUTO** setting has been activated in the Sipper options menu, the rinse agent is sucked in automatically once the measurement process is complete.

The zero measurement is shown on the display.

- Place the feed tube in the sample solution and press **READ**.


The sample solution is sucked in over the selected sip time. Once the settle time has passed, the measurement is taken.

**Note:** The remaining time is shown on the display. Press **CANCEL** to abort the operation.

**Note:** If the **PURGE START: AUTO** setting has been activated in the Sipper options menu, the rinse agent is sucked in automatically once the measurement process is complete.

**Note:** Either deionized water or the next sample can be used as a rinse agent.

The measurement is shown on the display.


## Execution of a Time Course Measurement

- Press **TIME COURSE**.
- Press **OPTIONS** to adjust the settings within time course mode.
- Place the feed tube in the zero solution and press **ZERO**.

The zero solution is sucked in over the selected sip time. Once the settle time has passed, the zero measurement is taken.

**Note:** If the **PURGE START: AUTO** setting has been activated in the Sipper options menu, the rinse agent is sucked in automatically once the measurement process is complete.

The zero measurement is shown on the display.

- Place the feed tube in the sample solution and press **READ**.  
The sample solution is sucked in over the selected sip time. Once the settle time has passed, the measurement is taken.  
The measurement is taken in accordance with the time course settings. The **ZERO/READ** buttons change to **MARK/STOP**. To abort the measurement process, press **STOP**.

Once the entire time has passed, the measurement is shown as a curve on the display.

**Note:** If the **PURGE START: AUTO** setting has been activated in the Sipper options menu, the rinse agent is sucked in automatically once the measurement process is complete.

## Maintenance

### Cleaning

# NOTICE

There may always be residual fluid in the tubes. Always wear suitable protective clothing and goggles to change the tubes or perform maintenance work in order to prevent contact with the skin or eyes.

### SIP 10 Module

Clean the SIP 10 module using a damp cloth. Do not use any type of solvent (for example, acetone).

### Pour-Thru Cell

Purge the Pour-Thru Cell before and after every series of tests with deionized water.

If the Pour-Thru Cell is heavily contaminated, repeat the purge cycle several times.

Before the Pour-Thru Cell is put into storage, perform the purge cycle several times with air to remove any residual fluid. Hold the Pour-Thru Cell with the OUT connection slanted downwards during this procedure.

Check the Pour-Thru Cell viewing window. If it is dirty or fogged up, place the Pour-Thru Cell in a soap solution or in dilute acid. Then rinse the Pour-Thru Cell thoroughly with deionized water.

### Tube Cleaning

Clean the tubes with deionized water after every series of measurements.


The tubes are exposed to chemicals and must be replaced periodically.

### Pump Tubing Replacement

The pump tubing is exposed to mechanical loads and chemicals and must be replaced periodically.

1. Perform the purge cycle several times with air to remove any residual fluid from the system.
2. Disconnect the USB cable from the photometer.
3. Remove the connection tube and discharge tube from the input and output connections on the peristaltic pump.
4. Tilt back the lever behind the peristaltic pump (Figure 2 step 1)  
The peristaltic pump cover will move back.
5. Lift up the peristaltic pump cover (Figure 2 step 2) and dispose of the cover, pump tubing and connections.
6. Position the new cover with pre-assembled pump tubing and connections on the peristaltic pump.
7. Push down the cover on the peristaltic pump and tilt up the lever behind the peristaltic pump (Figure 2 steps 3 and 4).  
The peristaltic pump cover will move forward.
8. Connect the connection tube and the discharge tube to the input and output connections on the peristaltic pump.
9. Connect the SIP 10 USB cable to a USB port on the photometer.

Figure 2 Pump tubing change

## Troubleshooting

### Error Messages

Potential SIP 10 errors are displayed by the photometer.

Table 3 Error Messages

Error displayed	Cause	Resolution
Sipper module not connected. Please check cable.	Connection between SIP 10 and photometer disrupted.	Check the USB cable. The cable length must not exceed 1 m (3.3 ft) and the cable must be directly connected to the photometer. Remove any other devices that are connected where necessary.
Please check sipper and tube.	Pump tubing not correctly inserted.	Loosen the cover and position again. Where required, run the pump briefly before the lever is moved to make sure that the tube is correctly positioned around the rollers.

## Replacement Parts

Designation	Catalog number
SIP 10 Sipper module set complete with tube set and 1 inch/1 cm Pour-Thru Cell, US	LQV157.99.10002
SIP 10 Sipper module; including pump tubing, US	LQV157.99.00002
USB cable, type AB 1 m (3.3 ft)	LZQ064
Dual pathlength Pour-Thru Cell, plastic, 1 inch/1 cm	LZV876
Complete tube set; including Tygon tubing, 1.70 m (5.6 ft), pump tubing (Lagoprene <sup>®</sup> ) pre-assembled with cover and connections, stainless steel tube adapter	LZV875
Pump tubing (Lagoprene) pre-assembled with cover and connections	LZV877
User manual, multilingual SIP 10, US	DOC012.97.90328
Information on Tygon tubes, EC no. 1907/2006 (REACH)	DOC063.98.90334

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© HACH Company, 2011. All rights reserved. Printed in Germany